

BIRTHING BED

VMS[®]
Virtue | Merit | Security



VBB-1001 Birthing Bed (Hydraulic Height Adjustment)

A **BOND** you can **TRUST**

BIRTHING BED



Product Specification:

Overall approx dimension: 2070mm (L) x 850mm (W) x 550mm to 850mm (H). Three section top. Backrest adjustable on screw mechanism. Leg section can be telescoped under the backrest for lithotomy positions. Removable stainless steel head & foot boards with laminated panels. Foot support for high/low chair positions. Middle & Leg section tops are removable for cleaning. A pair of upholstered mild steel lithotomy crutches mounted on SS rods. Complete with SS side railings, SS tray & SS telescopic IV pole. Height adjustable by foot operated hydraulic pump. Handle wound trendelenburg position. Three section mattress. Synthetic body castors of 125mm wheel dia, two with brake and two without brake. Pre-treated & powder coated finish. Supplied in semi knocked down condition.

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East),

Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customer@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.