



#### **VEX-1001 Examination Couch With Drawers & Cabinets**

Product Specification:

Overall approx dimension: 1890mm (L) x 560mm (W) x 840mm (H). Strong construction made from mild steel ERW rectangular tubes and machine pressed mild steel CRCA sheet. Fitted with rexine covered two sections cushioned top. Manually adjustable headrest section on gas spring mechanism. Upper section of the box construction with three separate drawers. Lower section of the box construction with three separate cabinets having doors with locking facility. Height adjustable B.P. apparatus tray provided near headrest. Built-in sliding step stool. Pre-treated and powder coated finish. Supplied in semi knocked down condition.

**VEX-1004 Examination-cum-Gynaec Couch With "U" cut, Lithotomy Rods, Plain Cushioned Top and Drawers & Cabinets**

Overall approx dimension:

1890mm (L) x 560mm (W) x 840mm (H). Strong construction made from mild steel ERW rectangular tubes and machine pressed mild steel CRCA sheet. Fitted with rexine covered plain cushioned top with removable "U" cut. Stainless steel bowl, a pair of adjustable stainless steel lithotomy rods with rexine straps. Upper section of the box construction with three separate drawers. Lower section of the box construction with three separate cabinets having doors with locking facility. Height adjustable B.P. apparatus tray provided near head side. Built-in sliding step stool. Pre-treated and powder coated finish. Supplied in semi knocked down condition.



**VEX-1005 Examination-Cum-Gynaec Table with Lithotomy Rods and with Cushioned Top**

Overall approx dimension:

1830mm (L) x 610mm (W) x 825mm (H). Strong mild steel ERW tubular frame work. Fitted with machine pressed double bent two section mild steel CRCA sheet top with removable "U" cut with fixed cushioned top. Manually adjustable headrest on ratchet mechanism. Stainless steel bowl, a pair of adjustable stainless steel lithotomy rods with rexine straps. Legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.



**VEX-1002 Examination Table With Cushioned Top**

Overall approx dimension:

1830mm (L) x 610mm (W) x 825mm (H). Strong mild steel ERW tubular frame work. Fitted with machine pressed double bent two section mild steel CRCA sheet top with fixed cushioned top. Manually adjustable headrest on ratchet mechanism. Legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.

\* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

\* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.

**Vinod Medical Systems Pvt. Ltd.**

**Corporate Office :** 523, Chintamani Plaza, Mohan Studio Compound, Chakala, Andheri – Kurla Road, Andheri (E), Mumbai – 400 099, India.

**Tel :** +91 22 42152020 / 42152021. **Fax :** +91 22 42152024 • **Email :** salesbfd@vinodmedical.com • **Web :** www.vinodmedical.com

**Manufacturing Plant:** C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171. Gujarat. India.