

FOOD TABLES



VBT-1001 Adjustable Bedside Table With Two Section Laminated Board Top

(Height adjustable on Gas Spring Mechanism)

Overall approx dimension: 810mm (L) x 355mm (W). Height adjustable from 760mm to 1050mm approx, with the help of gas spring. Fitted with two section laminated Board top with four side raised edges. Telescopic stem made from CRCA mild steel sheet. Top and base frame made from mild steel ERW tubes. Mounted on four synthetic body castors 50mm wheel dia, without brake. Pre-treated & powder coated finish. Supplied in knocked down condition.



VBT-1002 Adjustable Bedside Table With Two Section Laminated Board Top

(Height adjustable on Gear Mechanism)

Overall approx dimension: 810mm (L) x 355mm (W). Height adjustable from 760mm to 1050mm approx, with gear having stainless steel handle. Fitted with two section laminated Board top with four side raised edges. Telescopic stem made from CRCA mild steel sheet. Top and base frame made from mild steel ERW tubes. Mounted on four synthetic body castors 50mm wheel dia, without brake. Pre-treated & powder coated finish. Supplied in knocked down condition.



VBT-1004 Over Bed Table With Plain Laminated Board Top (Fixed Height)

Overall approx dimension: 1120mm (L) x 450mm (W) x 970mm (H). Top approx dimension: 760mm (L) x 450mm (W). Strong mild steel round tubular frame work. Mounted on four synthetic body castors 50mm wheel dia, without brake. Fitted with plain laminated board top. Pre-treated & powder coated finish. Supplied in knocked down condition.



VBT-1003 Adjustable Bedside Table With Plain SS Top (Height adjustable on Gear Mechanism)

Overall approx dimension: 810mm (L) x 355mm (W). Height adjustable from 760mm to 1050mm approx, with gear having stainless steel handle. Fitted with single section stainless steel table top with four side raised edges. Telescopic stem made from CRCA mild steel sheet. Top and base frame made from mild steel ERW tubes. Mounted on four synthetic body castors 50mm wheel dia, without brake. Pre-treated & powder coated finish. Supplied in knocked down condition.

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customer@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.