

PAEDIATRIC BEDS



VPB-1011 Paediatric Bed with Full Length Drop Type Side Railings

Overall approx dimension: 1480mm (L) x 840mm (W) x 600mm (H). Bed frame approx. dimension: 1370mm (L) x 760mm (W). Uniformly perforated mild steel CRCA sheet top. Both head bow and leg bow of equal height of approx 1060mm (H). A pair of full length drop type side railing. Legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.



VPB-1012 Baby Crib With Trolley

Overall approx dimension: 945mm (L) x 425mm (W) x 1030mm (H). Crib approx dimension: 760mm (L) x 380mm (W) x 300mm (H). Strong construction made of mild steel ERW round tubes and bright rods. Mounted on four synthetic body castors 125mm wheel dia, two with brake and two without brake. Pre-treated and powder coated finish. Supplied in knocked down condition.



VPB-1013 Baby Crib With Attachment

Crib approx dimension: 760mm (L) x 380mm (W) x 300mm (H). Crib construction made from mild steel ERW round tubes and bright rods. Swings on two vertical tubes to be fixed to leg end bow. Pre-treated and powder coated finish. Supplied in semi knocked down condition.

Optional Accessories for Paediatric Beds:

- VMT-1005 Rexine Covered Mattress Suitable For Paediatric Bed
- VMT-1006 Rexine Covered Mattress Suitable For Baby Crib

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customercare@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.