

RECOVERY BEDS



VRB-5002 Recovery Bed (Four Section)

(Removable SS head & foot boards with laminated panels)

Overall approx dimension: 2170mm (L) X 1000mm (W) X 600mm (H)



VRB-5001 Recovery Bed (Four Section)

(Removable Polymer moulded head & foot boards)

Overall approx dimension: 2240mm (L) x 1000mm (W) x 600mm (H)

RECOVERY BEDS (Four Section)

Product Specification:

Bed frame approx. dimension: 2080mm (L) x 910mm (W). Backrest, upper leg and Trendelenburg / Reverse Trendelenburg adjustments on crank mechanism with separate screws provided at leg end with stainless steel handles. Four section perforated mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Mounted on four synthetic body castors 125mm wheel dia, two with brake and two without brake. Telescopic stainless steel IV pole with four locations at four corners of the bed. Pre-treated and powder coated finish. Supplied in knocked down condition.

Optional Accessories For Recovery Beds:

- VAC-1002 Full Length Collapsible Side Railings With SS Tube & MS Frame (Pair)
- VAC-1003 Collapsible Side Railings With SS Tube & MS brackets (Pair)
- VAC-1004 Swing Away Type SS Tube Railings (Pair)
- VAC-1005 Swing Away Type MS Tube Railings (Pair)
- VAC-1007 Urine Bag Holder (MS EPC)
- VAC-1008 Chart Holder (MS EPC)
- VMT-1002 PU Foam Mattress – Two Sections
- VMT-1004 PU Foam Mattress – Four Sections

A **BOND** you can **TRUST**

RECOVERY BEDS



VRB-5004 Recovery Bed (Two Section)

(Removable SS head & foot boards with laminated panels)

Overall approx dimension: 2070mm(L) x 1000mm (W) X 600mm(H)



VRB-5003 Recovery Bed (Two Section)

(Removable Polymer moulded head & foot boards)

Overall approx dimension: 2140mm (L) x 1000mm (W) x 600mm (H)

RECOVERY BEDS (Two Section)

Product Specification:

Bed frame approx. dimension: 1980mm (L) x 910mm (W). Backrest and Trendelenburg/Reverse Trendelenburg adjustments on crank mechanism with separate screws provided at leg end with stainless steel handles. Two section perforated mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Mounted on four synthetic body castors 125mm wheel dia, two with brake and two without brake. Telescopic stainless steel IV pole with four locations at four corners of the bed. Pre-treated and powder coated finish. Supplied in knocked down condition.

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customercare@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.