

SHIFTING TROLLEY



VST-1002 Emergency & Recovery Trolley
(Height adjustable on foot operated indigenous hydraulic pump)

Product Specification:


Overall approx dimension: 1905mm (L) x 710mm (W). Stretcher approx dimension: 1830mm (L) x 555mm (W). Height adjustment approx: 660mm to 910mm. Height adjustable on foot operated indigenous hydraulic pump. Two section x-ray permeable laminated board top. Backrest adjustable on ratchet mechanism. Trendelenburg/Reverse Trendelenburg positions on hand operated gas spring mechanism. Mounted on four synthetic body castors 125mm wheel dia, two with brake and two without brake. Sliding X-ray cassette holder, utility tray, oxygen cylinder cage, corner rubber buffers & handles covered with synthetic material. A pair of swing away type stainless steel railings. Telescopic stainless steel IV pole with four locations. Pre-treated and powder coated finish. Supplied in semi knocked down condition.

Optional Accessories For Shifting Trolley

- VMT-1008 Two Section Rexine Covered Mattress Suitable For Emergency & Recovery Trolley
- VMT-1009 Single Section Rexine Covered Mattress Suitable For Stretcher Trolley
- VMT-1010 Rexine Covered Mattress Suitable For Folding Stretcher Trolley

A BOND you can TRUST

SHIFTING TROLLEY



VST-1004 Stretcher Trolley (MS Framework & SS Stretcher Top)

Overall approx dimension: 2060mm (L) x 560mm (W) x 820mm (H). Strong construction made from mild steel ERW round tubes. Mounted on four metallic body castors 200mm wheel dia, without brake. Removable stretcher top made of stainless steel sheet supported on mild steel tubular frame. Pre-treated and powder coated finish. Supplied in knocked down condition.



VST-1007 Stretcher Trolley - Folding (MS Framework & SS Stretcher Top)

Overall approx dimension:
2120mm (L) x 590 mm (W) x 820mm (H). Strong construction made from mild steel ERW round tubes. Mounted on four metallic body castors 200mm wheel dia, without brake. Fixed folding stretcher top made of stainless steel sheet supported on mild steel tubular frame. Pre-treated and powder coated finish. Supplied in knocked down condition.



VST-1005 Stretcher Trolley (SS Framework & SS Stretcher Top)

Overall approx dimension: 2060mm (L) x 560mm (W) x 820mm (H). Strong construction made from stainless steel round tubes. Mounted on four metallic body castors 200mm wheel dia, without brake. Removable stretcher top made of stainless steel sheet supported on stainless steel tubular frame. Supplied in knocked down condition.

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customercare@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.