

ORVEE



MEDICAL SOLUTIONS




BIS Certified

AERB Type Approved

Classic-100 Mobile X-Ray Radiography System





Advanced Technology:

Orvee Medical Solutions specializes in computer based diagnostic X-Ray Machines. The model Classic-100 has an attractive TFT monitor display which is embedded with an anatomical program of different body parts. This interface makes it appealing, simple and user friendly.

Mobility:

The perfectly balanced lightweight x-ray unit can be moved anywhere within the hospital at ease. The gas spring helps in positioning the tube to suit different bed heights. The stability is maintained by a well designed base with a balanced tube head.



REPRESENTED BY:

SPECIFICATIONS

- Generator** : 100mA and 100kV
- X-ray Tube** : Stationary Anode with a focal spot of 2.8mm or less.
- Tube Voltage** : 40kV to 100kV (Selection in steps of 1 kV)
- Tube Current** : 20mA, 40mA, 60mA, 80mA & 100mA
- Timer** : 0.02secs to 5secs solid state
- Maximum Power** : 8kW/10KVA with line resistance of 0.4 ohms or less

Manufactured By: **ORVEE MEDICAL SOLUTIONS PVT LTD, Bangalore-560027**

Marketed by :

Vinod Medical Systems Pvt. Ltd.

Certifications: ISO 9001:2008, ISO 14001:2004, ISO 13485:2003 and CE.

Regd. Off.: Kripa Kunj, B 1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, (C.G), India. Tel.: +91 771 4214400 / 21, Fax : +91 771 2882404

Corp. Off.: B - 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (E) Mumbai - 400 099. India. Tel.: +91 22 42152020 / 21, Fax.: +91 22 42152024

