

# ORVEE



MEDICAL SOLUTIONS



BIS Certified

AERB Type Approved

## Classic-100 R X-RAY RADIOGRAPHY SYSTEM





### Advanced Technology:

Orvee Medical Solutions specializes in computer based X-Ray Machines. The model Classic 100 R has an attractive TFT Monitor display which is embedded with an anatomical program of different body parts. This interface makes it Appealing, Simple & User Friendly.

### Salient Features

#### Anatomical Program

Enables users to select from more than 200 presets available

#### Radiographic Settings

Single step KVP & mAs settings for better Radiographic results

#### Automatic Line Voltage Compensation

### REPRESENTED BY:

#### GENERATOR

#### CLASSIC 100 R

Radiographic rating

100mA & 100KV

Tube Current

20mA, 40mA, 60mA, 80mA, 100mA

Radiographic KV

40 To 100KV, Variations In Steps Of 1KV

Exposure time

0.02 To 5.0 Seconds In Steps Of 0.02 Sec

X-ray tube

Single Focus Rotating Anode

Display

Computerized KV, mA, mAs & Time Indications

Overload Protection

Computer Controlled

LV compensation

15% Of Rated Input Voltage

Optimized Wiring Increases Efficiency Of Output Power

Manufactured By: **ORVEE MEDICAL SOLUTIONS PVT LTD**, Bangalore-560027

Marketed by :

### Vinod Medical Systems Pvt. Ltd.

Certifications: ISO 9001:2008, ISO 14001:2004, ISO 13485:2003 and CE.

**Regd. Off.:** Kripa Kunj, B 1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, (C.G), India. Tel.: +91 771 4214400 / 21, Fax : +91 771 2882404

**Corp. Off.:** B - 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (E) Mumbai - 400 099, India. Tel.: +91 22 42152020 / 21, Fax.: +91 22 42152024

