

ORVEE



MEDICAL SOLUTIONS




BIS Certified

AERB Type Approved

Classic-300 X-RAY RADIOGRAPHY SYSTEM





Advanced Technology:

Orvee Medical Solutions specializes in computer based X-Ray Machines. The model Classic 300 has an attractive TFT Monitor display which is embedded with an anatomical program of different body parts. This interface makes it Appealing, Simple & User Friendly.

Salient Features


Anatomical Program

Enables users to select from more than 200 presets available

Radiographic Settings

Single step KVP & mAs settings for better Radiographic results

Automatic Line Voltage Compensation



REPRESENTED BY:

GENERATOR	CLASSIC 300
Radiographic rating	300mA & 100KV
Small focus	50, 100mA
Large focus	150, 200, 250 & 300mA
Radiographic KV	40 To 100KV, Variations In Steps Of 1KV
Exposure time	0.02 To 5.0 Seconds In Steps Of 0.02 Sec
X-ray tube	Dual Focus Rotating Anode
Display	Computerized KV, mA, mAs & Time Indications
Overload Protection	Computer Controlled
LV compensation	15% Of Rated Input Voltage
Power requirement	3 Phase, 400/440 Volts, 30 KVa, 50 Cycles AC

Optimized Wiring Increases Efficiency Of Output Power

Manufactured By: **ORVEE MEDICAL SOLUTIONS PVT LTD, Bangalore-560027**

Marketed by :

Vinod Medical Systems Pvt. Ltd.

Certifications: ISO 9001:2008, ISO 14001:2004, ISO 13485:2003 and CE.

Regd. Off.: Kripa Kunj, B 1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, (C.G), India. Tel.: +91 771 4214400 / 21, Fax : +91 771 2882404

Corp. Off.: B - 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (E) Mumbai - 400 099. India. Tel.: +91 22 42152020 / 21, Fax.: +91 22 42152024

