

# ORVEE

MEDICAL SOLUTIONS






BIS Certified

AERB Type Approved

## Classic-500

### X-RAY RADIOGRAPHY SYSTEM





## Advanced Technology:

Orvee Medical Solutions specializes in computer based X-Ray Machines. The model Classic 500 has an attractive TFT Monitor display which is embedded with an anatomical program of different body parts. This interface makes it Appealing, Simple & User Friendly.

## Salient Features


### Anatomical Program

Enables users to select from more than 200 presets available

### Radiographic Settings

Single step KVP & mAs settings for better Radiographic results

### Automatic Line Voltage Compensation



**REPRESENTED BY:**

GENERATOR	CLASSIC 500
Radiographic rating	500mA & 100KV
Small focus	50, 100mA
Large focus	200, 300 & 500mA
Radiographic KV	40 To 100KV, Variations In Steps Of 1KV
Exposure time	0.02 To 5.0 Seconds In Steps Of 0.02 Sec
X-ray tube	Dual Focus Rotating Anode
Display	Computerized KV, mA, mAs & Time Indications
Overload Protection	Computer Controlled
LV compensation	15% Of Rated Input Voltage
Power requirement	3 Phase, 400/440 Volts, 50 KVa, 50 Cycles AC

Optimized Wiring Increases Efficiency Of Output Power

Manufactured By: **ORVEE MEDICAL SOLUTIONS PVT LTD, Bangalore-560027**

Marketed by :

## Vinod Medical Systems Pvt. Ltd.

Certifications: ISO 9001:2008, ISO 14001:2004, ISO 13485:2003 and CE.

**Regd. Off.:** Kripa Kunj, B 1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, (C.G), India. Tel.: +91 771 4214400 / 21, Fax : +91 771 2882404

**Corp. Off.:** B - 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (E) Mumbai - 400 099. India. Tel.: +91 22 42152020 / 21, Fax.: +91 22 42152024

