

# FOWLER BEDS




**VFB-4002 Fowler Bed**  
(Removable SS head & foot boards with laminated panels)  
2120mm (L) x 910mm (W) x 600mm (H)<sup>#</sup>



**VFB-4001 Fowler Bed**  
(Removable Polymer moulded head & foot boards)  
2240mm (L) x 955mm (W) x 600mm (H)<sup>#</sup>



**VFB-4003 Fowler Bed**  
(Removable MS head & foot boards )  
2120mm (L) x 910mm (W) x 600mm (H)<sup>#</sup>



**VFB-4004 Fowler Bed**  
(Un-equal MS head & foot boards)  
2180mm (L) x 910mm (W) x 600mm (H)<sup>#</sup>

## Fowler Beds

### Product Specification:

Bed frame approx dimension : 2080mm (L) x 910mm (W). Backrest, upper leg section adjusted on crank mechanism with separate screws provided at leg end with stainless steel handles. Calf section adjusted on ratchet mechanism. Uniformly perforated four section mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Four locations on bed frame for telescopic IV pole. "H" type legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.

### Optional Accessories for Fowler/Semi Fowler Beds:

- VAC-1001 Telescopic IV Pole SS
- VAC-1002 Full Length Collapsible Side Railings With SS Tube & MS Frame (Pair)
- VAC-1003 Collapsible Side Railings With SS Tube & MS brackets (Pair)
- VAC-1004 Swing Away Type SS Tube Railings (Pair)
- VAC-1005 Swing Away Type MS Tube Railings (Pair)
- VAC-1006 Mattress Arrestor ( set of four)
- VAC-1007 Urine Bag Holder (MS EPC)
- VAC-1008 Chart Holder (MS EPC)
- VAC-1010 Castor wheel ( a set of four , 125mm wheel dia synthetic body castors , two with brake and two without brake )
- VMT-1002 PU Foam Mattress – Two Sections
- VMT-1004 PU Foam Mattress – Four Sections