

ICU BEDS



VIB-6002 Intensive Care Bed (Manual)

(Removable SS head & foot boards with laminated panels)

Overall approx dimension:

2170mm (L) x 1000 mm (W)

Product Specification:

Bed frame approx. dimension: 2080mm (L) x 910mm (W). Height adjustment approx. (without mattress): Minimum height 490mm to Maximum height 710mm. Backrest, upper leg, height and Trendelenburg / Reverse Trendelenburg adjustments on crank mechanism with separate screw mechanisms provided at leg end with stainless steel handles. Four section perforated mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Mounted on four synthetic body castors 125mm wheel dia, two with brake and two without brake. Rubber buffers at four corners. Telescopic stainless steel IV pole with four locations at four corners of the bed. Pre-treated and powder coated finish. Supplied in knocked down condition.

A **BOND** you can **TRUST**

ICU BEDS



VIB-6001 Intensive Care Bed (Manual)

(Removable Polymer moulded head & foot boards)

Overall approx dimension:

2240mm (L) x 1000 mm (W)

Optional Accessories for ICU Beds:

- VAC-1002 Full Length Collapsible Side Railings With SS Tube & MS Frame (Pair)
- VAC-1003 Collapsible Side Railings With SS Tube & MS brackets (Pair)
- VAC-1004 Swing Away Type SS Tube Railings (Pair)
- VAC-1005 Swing Away Type MS Tube Railings (Pair)
- VAC-1007 Urine Bag Holder (MS EPC)
- VAC-1008 Chart Holder (MS EPC)
- VMT-1004 PU Foam Mattress - Four Sections

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customercare@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.