

OBSTETRICS

VMS[®]
Virtue | Merit | Security



VLT-1001 Obstetric Labour Table
(SS Framework & SS Table Top with Lithotomy Rods)

Product Specification:

Overall approx dimension: 1880mm (L) x 900mm (W) x 800mm (H) (fully extended). Stainless steel table top in three sections. Stainless steel tubular framework. Backrest adjustable on ratchet mechanism. Middle section with "U" cut. Leg section can be slid underneath main frame for lithotomy positions. Trendelenburg position on hand operated gas spring mechanism. Removable stainless steel head and side railings, stainless steel hand grips and stainless steel bowl. A pair of rubber padded stainless steel lithotomy rods with rexine straps. Supplied in semi knocked down condition.

A BOND you can TRUST



VLT-1003 Obstetric Labour Table
(MS Framework & SS Table Top with Lithotomy Rods)

Optional Accessories for Obstetric Table:

- VAC-1001 Telescopic IV Pole SS
- VMT-1007 Rexine Covered Mattress Suitable For Obstetric Labour Table
- VAC-1013 Padded Lithotomy Crutches Suitable For Obstetric Labour Table



Lithotomy Crutches

Product Specification:

Overall approx dimension: 1880mm (L) x 900mm (W) x 800mm (H) (fully extended) . Stainless steel table top in three sections . Mild steel tubular framework . Backrest adjustable on ratchet mechanism. Middle section with "U" cut. Leg section can be slid underneath main frame for lithotomy positions. Trendelenburg position on hand operated gas spring mechanism. Removable stainless steel head and side railings, stainless steel hand grips and stainless steel bowl. A pair of rubber padded stainless steel lithotomy rods with rexine straps. Supplied in semi knocked down condition. Pre-treated and powder coated finish.

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.
* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.

Vinod Medical Systems Pvt. Ltd.

Corporate Office : 523, Chintamani Plaza, Mohan Studio Compound, Chakala, Andheri – Kurla Road, Andheri (E), Mumbai – 400 099, India.
Tel : +91 22 42152020 / 42152021. **Fax :** +91 22 42152024 • **Email :** salesbfd@vinodmedical.com • **Web :** www.vinodmedical.com
Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate , Near Apsara Pencil, Umbergaon - 396 171. Gujarat. India.