

SEMI FOWLER BEDS



VSB-3002 Semi Fowler Bed
(Removable SS head & foot boards
with laminated panels)
2020mm (L) x 910mm (W) x 600mm (H)[#]



VSB-3001 Semi Fowler Bed
(Removable Polymer moulded head & foot boards)
2140mm (L) x 955mm (W) x 600mm (H)[#]



VSB-3003 Semi Fowler Bed
(Removable MS head & foot boards)
2020mm (L) x 910mm (W) x 600mm (H)[#]



VSB-3004 Semi Fowler Bed
(Un-equal MS head & foot boards)
2080mm (L) x 910mm (W) x 600mm (H)[#]

Semi Fowler Beds Product Specification:

Bed frame approx dimension: 1980mm (L) x 910mm (W). Backrest adjustable on crank mechanism with separate screw provided at leg end with stainless steel handle. Uniformly perforated two section mild steel CRCA sheet top. Framework made from mild steel ERW rectangular tubes. Four locations on bed frame for telescopic IV pole. "H" type legs fitted with rubber shoes. Pre-treated and powder coated finish. Supplied in knocked down condition.

Vinod Medical Systems Pvt. Ltd.

Regd. Office: Kripa Kunj, B1-2, Sai Nagar, Near Railway Crossing, Raipur - 492 009, Chhattisgarh. Tel.: +91 771 4214400 / 21 / 22

Corporate Office: B / 523, Chintamani Plaza, Mohan Studio Compound, Andheri Kurla Road, Chakala, Andheri (East), Mumbai - 400 099, Maharashtra. Tel.: +91 22 42152020 / 21, Fax: +91 22 42152024, Email: customercare@vinodmedical.com

Manufacturing Plant: C-1-275 / 16 & 17, G.I.D.C. Estate, Near Apsara Pencil, Umbergaon - 396 171, Gujarat.

Website: www.vinodmedical.com

* Accessories and colours shown in the catalogue illustrations may not be the part of standard product configuration.

* We reserve the right to make changes without prior notice in design, specifications and illustrations. Reproduction of this product catalogue in any form is liable for legal action.